

TalkTools Jiggler® Combo - Dragon & Frog



Price: €48.80 inc. 19% VAT

SKU: 103025

Product Categories: [ORAL MOTOR](#), [VIBRATING TOOLS & ACCESSORIES](#)

Product Page:

<https://theraplaycy.com/product/talktools-jiggler-combo-dragon-frog/>

Product Description

All four TalkTools Jiggler® are effective for increasing sensory awareness in the muscles of the jaw, lips, and reducing "fixing" or tightening in the muscles of the jaw. In some cases the Jigglers, satisfy the need for temporomandibular joint (TMJ) stimulation as a replacement to negative oral habits. In addition, each individual TalkTools Jiggler® can be used to meet more specific needs:

Dragon Therapeutic Goals:

GOAL #1: To increase sensory awareness.

GOAL #2: To improve lip rounding and tongue retraction.

Feeding - Straw Drinking

OPT - Bubble blowing and round mouth horns

Speech - /oo/ /oh/ /w/

GOAL #3: To improve lip rounding and lip protrusion with the added texture on the Dragon for clients who need additional sensory input.

Feeding -

Straw Drinking

OPT - Bubble blowing and round mouth horns

Speech - /oo/ /oh/ /w/ /r/ /sh/ /ch/

Frog Therapeutic Goals:

GOAL #1: To reduce "fixing" (tightening in the muscles of the jaw) before muscle-based exercises, feeding and speech. This will satisfy the need for TMJ stimulation and release of a tonic bite.

GOAL #2: To increase sensory awareness, using a smooth surface, in the muscles of the jaw and lips prior to the introduction of muscle-based exercises.

Feeding - Spoon feeding and cup drinking

OPT - Flat mouth horns, tongue depressor activity for lip closure

Speech - /m/ /b/ /p/

GOAL #3: To improve lip rounding and lip protrusion.

Feeding - Straw Drinking

OPT - Bubble blowing and round mouth horns

Speech - /oo/ /oh/ /w/ /r/ /sh/ /ch/

CONTRADICTIONS: Do NOT use this tool with clients who have been diagnosed with seizure disorders that are stimulated by vibration.

☐☐ These tools are intended for external use with adult supervision. Chewing on these tools is not recommended.

Requires 2 AA Batteries which are NOT included.

Use this tool for /oo/ /oh/ /w/ /r/ /sh/ /ch/ /m/ /b/ /p/ sounds

DETAILS

test