

TalkTools Jiggler® Combo - Unicorn & Sloth



Price: €48.80 inc. 19% VAT

SKU: 103024

Product Categories: [ORAL MOTOR](#), [VIBRATING TOOLS & ACCESSORIES](#)

Product Page:

<https://theraplaycy.com/product/talktools-jiggler-combo-unicorn-sloth/>

Product Description

All four TalkTools Jiggler® are effective for increasing sensory awareness in the muscles of the jaw, lips, and reducing "fixing" or tightening in the muscles of the jaw. In some cases the Jigglers, satisfy the need for temporomandibular joint (TMJ) stimulation as a replacement to negative oral habits. In addition, each individual TalkTools Jiggler® can be used to meet more specific needs:

TalkTools Jiggler® Combo - Unicorn & Sloth

Unicorn Therapeutic Goals:

GOAL #1: To reduce "fixing" (tightening in the muscles of the jaw) before muscle-based exercises, feeding and speech. This will satisfy the need for TMJ stimulation and release of a tonic bite.

GOAL #2: To improve lip closure and tongue retraction before muscle-based exercises, feeding and speech.

Feeding - Spoon feeding and cup drinking

OPT - Flat mouthed horns & tongue depressor

Speech - /m/ /b/ /p/

GOAL #3: To introduce spoon feeding for clients including those with oral aversions.

Sloth Therapeutic Goals:

GOAL #1: To reduce "fixing" (tightening in the muscles of the jaw) prior to the introduction of muscle-based exercises and/or to satisfy the need for TMJ stimulation.

GOAL #2: To improve sensory awareness in the muscles of the jaw and lips, prior to the introduction of muscle-based exercises.

Feeding - Spoon Feeding

OPT - Flat mouthed horns and tongue depressor

Speech - /m/ /b/ /p/

GOAL #3: To improve lip rounding.

Feeding - Straw drinking

OPT - Bubble blowing and round mouth horns

Speech - /oo/ /oh/ /w/

CONTRADICTIONS: Do NOT use this tool with clients who have been diagnosed with seizure disorders that are stimulated by vibration.

☐☐ These tools are intended for external use with adult supervision. Chewing on these tools is not recommended.

Requires 2 AA Batteries which are NOT included.

Use this tool for /oo/ /oh/ /w/ /r/ /sh/ /ch/ /m/ /b/ /p/ sounds

DETAILS

test